

MediNavix™

Guidewire





MediNavix™ Guide Wire

Precision Navigation. Reliable Performance.

MediNavix™ Guide Wire is designed for smooth navigation through coronary vasculature, enabling precise placement of angiography catheters, balloon catheters, and guiding catheters during interventional procedures. The hydrophilic-coated surface enhances trackability while providing excellent torque response and procedural control.

Key Features

- » Distal side Hydrophilic coating for reduced friction and smooth lesion crossing.
- » Excellent torque transmission and tactile feedback.
- » Radiopaque platinum-nickel spring tip for enhanced visibility.
- » Distal part Nitinol core technology for superior flexibility and shape retention.
- » Proximal part PTFE-coated stainless-steel core for optimal support and pushability.
- » Available in Straight and J-Tip (45°) configurations.
- » EO sterilized and single-use design.

Product Specifications

- » Guidewire O.D. : 0.014" (0.36 mm)
- » Guidewire Lengths: 180 cm, 190 cm and 300 cm
- » Tip Shapes: Straight / J-Tip
- » Coating: Hydrophilic Tip and PTFE Shaft Design

Guidewire length (cm)	Platinum-nickel spring ring Length (cm)	Total Spring ring length (cm)	Guidewire Tip shape	Catalogue number
180	2.5	2.5	straight	MNNSI4A25025S
180	2.5	2.5	J	MNNSI4A25025J
180	3	3	straight	MNNSI4A03030S
180	3	3	J	MNNSI4A03030J
180	3	11	straight	MNNSI4A03110S
180	3	11	J	MNNSI4A03110J
180	3	12	straight	MNNSI4A03120S
180	3	12	J	MNNSI4A03120J
180	3	20	straight	MNNSI4A03200S
180	3	20	J	MNNSI4A03200J
180	3	28	straight	MNNSI4A03280S
180	3	28	J	MNNSI4A03280J
190	3	4.5	straight	MNNSI4B03045S
190	3	4.5	J	MNNSI4B03045J
300	2.5	2.5	straight	MNNSI4C25025S
300	2.5	2.5	J	MNNSI4C25025J
300	3	3	straight	MNNSI4C03030S
300	3	3	J	MNNSI4C03030J
300	3	4.5	straight	MNNSI4C03045S
300	3	4.5	J	MNNSI4C03045J
300	3	11	straight	MNNSI4C03110S
300	3	11	J	MNNSI4C03110J
300	3	12	straight	MNNSI4C03120S
300	3	12	J	MNNSI4C03120J
300	3	20	straight	MNNSI4C03200S
300	3	20	J	MNNSI4C03200J
300	3	28	straight	MNNSI4C03280S
300	3	28	J	MNNSI4C03280J



MedivationBio

Achieving Possibilities

Medivation Bio is an emerging medical device company committed to delivering safe, reliable and high-quality consumables for critical care and routine medical use. Our product portfolio includes IV Cannulas, Hypodermic Syringes, Blood Collection Tubes, Chest Drainage Catheters, High-Pressure Tubings and Cardiovascular Devices.

We follow strict quality protocols and design products on:

Patient safety

Reliable engineering standards

Ease of use

Clinical performance

"Achieving Possibilities"

is our commitment to creating devices that support clinicians and improve patient care.

SIDCO Industrial Complex, Bari Brahmana EPIP Industrial Estate,
Samba, Jammu and Kashmir, India 181133.